

OTS15 Spacious Sitting Chamber

Model OTS15 · specification sheet

Spacious upright sitting chamber at 1.5 ATA with panoramic windows and room to read, work or relax.



Model	OTS15
Build type	Soft shell
Working pressure	1.5 ATA · 7 psi
Dimensions	Spacious upright, folding chair
Capacity	Sit, 1 person
Warranty	1 year (extendable +2 years)
Lead time	10–30 days (after testing)
Shipping	Sea / air freight, worldwide
MOQ / OEM	1 unit · OEM, ODM & colour customise

Key features

- Spacious upright interior
- 1.5 ATA (7 PSI) comfort
- Large panoramic windows
- Included folding chair

Overview

The Oxygen Times OTS15 is designed for people who appreciate comfort, space, and premium craftsmanship. Featuring a spacious upright design and operating at 1.5 ATA, it creates a relaxing environment where wellness becomes part of your daily lifestyle. Whether you're reading, working on a laptop, watching your favorite show, or simply enjoying quiet time, the OTS15 lets you experience hyperbaric oxygen sessions in exceptional comfort. Perfect for home wellness spaces, fitness studios, wellness centers, beauty salons, luxury resorts, and personal wellness rooms.

Technical specifications

Chamber Type	Sitting Hyperbaric Oxygen Chamber
Operating Pressure	Up to 1.5 ATA (7 PSI)
Chamber Material	High-Strength TPU + Protective Nylon Cover
Oxygen Concentration	Up to 93% ±3%
Oxygen Delivery	Mask or Nasal Cannula
Sealing System	Reinforced Triple-Zipper

Air Supply	Quiet Oil-Free Compressor
Safety Features	Dual Emergency Release Valves, Pressure Gauge

What ships in the box

- OTS15 Chamber
- Integrated Air & Oxygen System
- Oxygen Mask
- Nasal Cannula
- Air Hoses
- Pressure Gauge
- Repair Kit
- User Manual
- Power Cable

Quality control before shipment

- Pressure testing before dispatch
- Airtightness inspection
- Functional performance testing
- Appearance inspection
- Export packaging inspection

Certifications

TBC — no certificate numbers or issuing bodies are published for this model. Request the current certification position for your market from the factory before ordering.

Ordering

Minimum order	1 unit, including colour customisation
Lead time	10-30 days after quality testing, confirmed at order
Freight	Air freight 7-10 days, sea freight 4-6 weeks
OEM / ODM	Logo, colour, packaging, manuals, accessories, pressure configuration, destination electrical standard
Warranty	1 year, extendable by 2 further years
Contact	sales@o2times.com · WhatsApp +86 180 6880 3928

This document describes equipment. It is not medical advice. How a hyperbaric chamber may be used, and by whom, is set by the regulations of the country it is installed in, and meeting them is the buyer's responsibility. No product certification is claimed; values the factory does not publish are marked TBC.