

OTS20 2.0 ATA Soft Sitting Chamber

Model OTS20 · specification sheet

Professional portable sitting chamber delivering a true 2.0 ATA in high-strength TPU with a foldable design.



Model	OTS20
Build type	Soft shell
Working pressure	2.0 ATA
Dimensions	Foldable TPU sitting chamber
Capacity	Sit, 1 person
Warranty	1 year (extendable +2 years)
Lead time	10–30 days (after testing)
Shipping	Sea / air freight, worldwide
MOQ / OEM	1 unit · OEM, ODM & colour customise

Key features

- True 2.0 ATA performance
- Professional-grade TPU + YKK airtight zipper
- Reinforced internal support frame
- Foldable for transport and storage

Overview

2.0 ATA Soft Sitting Chamber The Oxygen Times OTS20 is a professional portable hyperbaric chamber engineered to deliver stable pressures of up to 2.0 ATA. Built with high-strength TPU, a reinforced support frame, and multiple integrated safety systems, it combines professional performance with the convenience of a foldable design. Designed for users who require true 2.0 ATA mild hyperbaric oxygen therapy, the OTS20 is ideal for wellness centers, rehabilitation facilities, sports recovery, and premium home installations. Request a Quote Email About OTS20 Oxygen Times provides complete OEM and ODM manufacturing services. Available customization includes: Company logo and branding Chamber and protective cover colors Operating pressure configuration Control system Oxygen supply configuration Accessories package Packaging design User manuals in multiple languages

What ships in the box

Every OTS20 package includes the essential equipment required for operation, eliminating the need for additional purchases. Dimensions: 220 × 70 × 110 cm The spacious chamber allows users to sit upright comfortably while also providing enough room to recline or lie down. A reinforced internal support frame ensures excellent structural stability during operation. An external support frame is also available upon request. The all-in-one system combines an air compressor and oxygen concentrator into a single compact unit, delivering oxygen concentrations of up to 93% for consistent and reliable performance. Includes multiple oxygen masks and nasal cannulas to accommodate different user preferences and operating requirements. Supplied with all required hoses, valves, connectors, fittings,

maintenance accessories, and comprehensive operating manuals for quick installation. Our experienced technical team provides remote installation guidance, commissioning assistance, operation training, and after-sales technical support to ensure a smooth setup.

Quality control before shipment

- Pressure testing before dispatch
- Airtightness inspection
- Functional performance testing
- Appearance inspection
- Export packaging inspection

Certifications

TBC — no certificate numbers or issuing bodies are published for this model. Request the current certification position for your market from the factory before ordering.

Ordering

Minimum order	1 unit, including colour customisation
Lead time	10-30 days after quality testing, confirmed at order
Freight	Air freight 7-10 days, sea freight 4-6 weeks
OEM / ODM	Logo, colour, packaging, manuals, accessories, pressure configuration, destination electrical standard
Warranty	1 year, extendable by 2 further years
Contact	sales@o2times.com · WhatsApp +86 180 6880 3928

This document describes equipment. It is not medical advice. How a hyperbaric chamber may be used, and by whom, is set by the regulations of the country it is installed in, and meeting them is the buyer's responsibility. No product certification is claimed; values the factory does not publish are marked TBC.