

OTSA15 Inflatable Lying Chamber

Model OTSA15 · specification sheet

Compact inflatable lying chamber with an air-column support system — lightweight, quick to set up, stable at 1.5 ATA.



Model	OTSA15
Build type	Soft shell
Working pressure	1.5 ATA
Dimensions	Air-column supported, Ø80 × 200 cm class
Capacity	Lie down, 1 person
Warranty	1 year (extendable +2 years)
Lead time	10–30 days (after testing)
Shipping	Sea / air freight, worldwide
MOQ / OEM	1 unit · OEM, ODM & colour customise

Key features

- Reinforced inflatable air-column frame
- Stable 1.5 ATA operation
- Lightweight and portable
- Quick setup, convenient storage

Overview

The Oxygen Times OTSA15 is a compact inflatable soft-shell hyperbaric chamber designed to deliver a comfortable lie-flat experience with reliable 1.5 ATA operating pressure. Unlike conventional soft chambers, the OTSA15 features an innovative inflatable air-column support system that enhances structural stability while maintaining excellent portability. Lightweight, easy to install, and engineered for dependable daily operation, it is an ideal solution for home wellness spaces, recovery studios, fitness facilities, and mobile applications. Since its introduction, the OTSA15 platform has maintained an outstanding record for reliable performance, reflecting Oxygen Times' commitment to quality and safety.

Technical specifications

Chamber Type	Inflatable Soft-Shell Lying Chamber
Operating Pressure	Up to 1.5 ATA (7 PSI)
Chamber Material	High-strength TPU with reinforced nylon exterior
Structure	Inflatable Air-Column Support
Sealing System	Reinforced Airtight Zipper

Oxygen Supply	Up to 93% via Mask or Nasal Cannula
Lead Time	10–30 Days After Quality Testing
OEM / ODM	Available

Quality control before shipment

- Pressure testing
- Airtightness inspection
- Functional performance testing
- Appearance inspection
- Secure export packaging for air or sea freight

Certifications

TBC — no certificate numbers or issuing bodies are published for this model. Request the current certification position for your market from the factory before ordering.

Ordering

Minimum order	1 unit, including colour customisation
Lead time	10-30 days after quality testing, confirmed at order
Freight	Air freight 7-10 days, sea freight 4-6 weeks
OEM / ODM	Logo, colour, packaging, manuals, accessories, pressure configuration, destination electrical standard
Warranty	1 year, extendable by 2 further years
Contact	sales@o2times.com · WhatsApp +86 180 6880 3928

This document describes equipment. It is not medical advice. How a hyperbaric chamber may be used, and by whom, is set by the regulations of the country it is installed in, and meeting them is the buyer's responsibility. No product certification is claimed; values the factory does not publish are marked TBC.