

OTSE15 Egg-Shaped Sitting Chamber

Model OTSE15 · specification sheet

Extra-spacious egg-shaped chamber that lets one adult sit upright, recline or lie down at a stable 1.5 ATA.



Model	OTSE15
Build type	Soft shell
Working pressure	1.5 ATA
Dimensions	Egg-shaped, extra headroom
Capacity	Sit, recline or lie
Warranty	1 year (extendable +2 years)
Lead time	10–30 days (after testing)
Shipping	Sea / air freight, worldwide
MOQ / OEM	1 unit · OEM, ODM & colour customise

Key features

- Egg-shaped extra-headroom interior
- Sit, recline or lie freely
- Internal shape-retaining frame
- Stable 1.5 ATA operation

Overview

The Oxygen Times OTSE15 is designed for users who value comfort, space, and ease of use. Unlike conventional horizontal chambers, the OTSE15 features an innovative egg-shaped structure with a spacious interior that allows users to sit upright, relax, recline, or lie down throughout the session. Built for daily wellness routines, premium wellness spaces, and residential environments, the OTSE15 combines a stable 1.5 ATA operating pressure with advanced safety engineering, durable construction, and effortless operation.

Technical specifications

Chamber Type	Extra-Spacious Sitting Chamber
Operating Pressure	Up to 1.5 ATA (7 PSI)
Chamber Material	High-strength TPU with reinforced nylon exterior
Structure	Internal Support Frame
Sealing System	Triple Airtight Zipper
Viewing Windows	Large Transparent Windows

Oxygen Supply	Up to 93% via Mask or Nasal Cannula
Included Seating	Adjustable Folding Chair
Safety Features	Manual Valve, Automatic Pressure Relief, Internal Pressure Gauge

Quality control before shipment

- Pressure testing before dispatch
- Airtightness inspection
- Functional performance testing
- Appearance inspection
- Export packaging inspection

Certifications

TBC — no certificate numbers or issuing bodies are published for this model. Request the current certification position for your market from the factory before ordering.

Ordering

Minimum order	1 unit, including colour customisation
Lead time	10-30 days after quality testing, confirmed at order
Freight	Air freight 7-10 days, sea freight 4-6 weeks
OEM / ODM	Logo, colour, packaging, manuals, accessories, pressure configuration, destination electrical standard
Warranty	1 year, extendable by 2 further years
Contact	sales@o2times.com · WhatsApp +86 180 6880 3928

This document describes equipment. It is not medical advice. How a hyperbaric chamber may be used, and by whom, is set by the regulations of the country it is installed in, and meeting them is the buyer's responsibility. No product certification is claimed; values the factory does not publish are marked TBC.